

# Sectra IDS7 BJC Healthcare

## Shortcuts for Radiologists

### Mouse Configuration

| Mouse button               | Image  | Stack             |
|----------------------------|--|-------------------|
| Left                       | Zoom   | Browse            |
| Right                      | Right-click toolbox  |                   |
| Middle                     | Window/Level   | Window/Level      |
| Left + Middle              | Zoom   | Zoom              |
| Left + Right               | Pan  | Pan               |
| Middle + Right             | Free Rotation  | Free Rotation     |
| Mouse wheel                | <None>   | Browse            |
| Right + scroll mouse wheel | Left/Right Series  | Left/Right Series |
| Middle button double click | Cine loop  |                   |
| Left Double-click          | Zoom one partition to full screen  |                   |
| Ctrl Alt Left-click        | Move series between partitions. Thumbnail will be highlighted = know which series is in the layout |                   |
| Ctrl + Scroll mouse        | Zoom in/out to cursor  |                   |
| Shift + Double Click       | Set Crosshair  |                   |

### General shortcuts image window

| Command                 | Shortcut |
|-------------------------|----------|
| Show Thumbnails         | T        |
| Show document (Reports) | R        |
| Reference Guide         | F1       |
| Refresh worklist        | F5       |

### Maneuvering between series/layouts

| Command              | Shortcut |
|----------------------|----------|
| Display Left Series  | ←        |
| Display Right Series | →        |
| Next Layout          | NUM .    |
| Previous Layout      | NUM 0    |
| Park                 | P        |

### Manipulating images

| Command                      | Shortcut     |
|------------------------------|--------------|
| Text                         | N            |
| Distance                     | D            |
| Orthogonal Distance          | O            |
| ROI circle                   | H            |
| Angle                        | Alt A        |
| Magnifying Glass             | M            |
| Zoom All Images to this size | Ctrl+Shift+Z |
| Invert                       | I            |
| Sharpen Edges                | Ctrl E       |
| Flip Horizontal              | Ctrl H       |
| Flip Vertical                | Ctrl V       |
| Reset                        | Ctrl+R       |
| Rotate 90° CW                | Alt+R        |
| Sync Group 5                 | NUM +        |
| No Sync Group                | NUM -        |
| Arrow                        | A            |

### Window/level presets

| Command           | Shortcut |
|-------------------|----------|
| Next window level | Ctrl+W   |
| Bone              | 1        |
| Head              | 2        |
| Stroke            | 3        |
| Liver             | 4        |
| Lung              | 5        |

|                           |   |
|---------------------------|---|
| Angio                     | 6 |
| Abdomen                   | 7 |
| Min/Max                   | 8 |
| Show window/Level Control | W |

### Stack handling

| Command   | Shortcut            |
|---|---------------------|
| Show/hide Crosshair                                 | X                   |
| or: Shift + Double-click                            |                     |
| Show/hide Localizer                                 | L                   |
| Show/hide Scanogram                                 | Shift+T             |
| First slice   | Home                |
| Last slice  | End                 |
| Switch slice by slice (shows each slice in a stack) | Page Down / Page Up |
| Scroll between series                               | Alt+S               |
| Automatically play cine loop                        | Ctrl L              |
| or: Middle button double-click                      |                     |
| Show Cine Controller                                | C                   |
| Stop All Cine Loops                                 | Ctrl +Shift+ K      |
| Set Key Image (single slice)                        | NUM /               |
| Reset Key Image                                     | NUM *               |
| Go to Previous Key Image                            | NUM 9               |
| Go to Next Key Image                                | NUM 3               |
| Slab  | -                   |
| Slab Properties                                     | Ctrl -              |
| MPR   | Shift+M             |
| Axial MPR   | Shift+A             |
| Coronal MPR   | Shift+C             |
| Sagittal MPR  | Shift+S             |
| 3D  | Shift+3             |
| Close 3D/MPR mode                                   | Esc                 |



